Handleiding

voor

Wenoll-System WS 100

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




3 APPLICATION


4 O₂ DOSAGE


One patient with Wenoll re-circulation system

5 MANUAL VENTILATION

6 APPLICATION WITH DEMAND MODULE (Option)

WS 100 Emergency Oxygen Case

The WS 100 is Wenoll's most compact emergency oxygen case of the hard cover - case series. It contains a 2-litre O₂ cylinder filled with medical oxygen and the CS 215 pressure regulator and is delivered assembled. It permits additional applications and therapy plans, e.g. one patient for **7.5 hours** or two patients simultaneously for 3 hours with **100% oxygen** (closed re-circulation system). The CS 215 pressure regulator can be adjusted in steps and has two additional medium pressure outlets (click connection) for the connection of **demand systems** (optional). The O₂ mask with reservoir bellows (open system) or the Wenoll re-circulation system is connected to the universal tubing nipple. The breathing system is cased ready to use and can be applied to the patient via a breathing mask or a mouthpiece.

The case is made from high-impact pressure resistant ABS plastic (polycarbonate) with a hermetically sealing O-ring. It is 100% water and dust tight, absolutely suitable for tropical use, is extremely stable under pressure and has air pressure compensation by means of a bleed system built into the locking arrangement. The upper part of the specially developed case inserts offers additional free space for individual stocks of emergency equipment.

■ Contents WS 100 with CS 215 Pressure Regulator

- | | |
|---|--|
| 1 off, case, yellow with pictogram and diving accident management | 1 off, breathing mask with air cushion seal |
| 2 off, case inserts, white | 1 off, mask holder for breathing mask |
| 1 off, O ₂ cylinder, 2 litre, filled with medical O ₂ , 200 bar | 1 off, mouthpiece |
| 1 off, CS 215 pressure regulator, 0-15 litres/min. in steps | 1 off, nose clip |
| 1 off, CO ₂ absorber with O ₂ connecting tube* | 1 off, O ₂ mask with reservoir bellows |
| 1 off, tubing set for WS 100/200/300 breathing system* | 1 off, APL valve, 0-60 mbar for manual respiration |
| 1 off, respiratory return bellows, 3 litres* | 1 off, one-way syringe for breathing mask |
| | 1 off, user manual and various product brochures |

■ Contents WS 100 with CS 215 Pressure Regulator and Demand Module

- | | |
|---|--|
| 1 off, case, yellow with pictogram and diving accident management | 1 off, breathing mask with air cushion seal |
| 2 off, case inserts, white | 1 off, mask holder for breathing mask |
| 1 off, O ₂ cylinder, 2 litre, filled with medical O ₂ , 200 bar | 1 off, mouthpiece |
| 1 off, CS 215 pressure regulator, 0-15 litres/min in steps | 1 off, nose clip |
| 1 off, demand module including connecting tube | 1 off, O ₂ mask with reservoir bellows |
| 1 off, CO ₂ absorber with O ₂ connection tube* | 1 off, APL valve, 0-60 mbar for manual respiration |
| 1 off, tubing set for WS 100/200/300 breathing system* | 1 off, one-way syringe for breathing mask |
| 1 off, respiratory return bellows, 3 litres* | 1 off, user manual and various product brochures |

■ WS 100 - Dimensions and Weights

Internal dimensions: 450 x 325 x 130 mm (BxDxH)

Internal volume: 19 litre

Weight, complete: 9.8 kg

■ Ordering Information

WS 100 with CS 215 pressure regulator _____ Art.-Nr. 121 015

WS 100 with CS 215 pressure regulator and demand module _____ Art.-Nr. 121 016

* Wenoll re-circulation system, complete, ready to run

Pressure Regulator CS 215 for med. O₂ / 200 bar

CE
0482

The CS 215 pressure regulator has been designed, manufactured and certified according to the newest guidelines DIN EN ISO 10524-1 (CE 0482). Its range of adjustment and functional components have been chosen such that all current O₂ applications can be used:

- constant flow systems
(eg. nasal cannula)
- open systems
(O₂ mask with or without reservoir bellows)
- demand systems
(pneumatic or electronic)
- closed systems
(re-circulation systems)

The patented flush function (40 litres/min.) with automatic reset permits fast flushing or filling of the breathing system for high O₂ concentrations.

The CS 215 pressure regulator has two self-sealing medium pressure outlets (click connection) and a constant-flow outlet with a universal tubing nipple for connecting tubes with a diameter of 4-6 mm.

The use of lightweight construction technology has given an additional weight reduction. Compared to conventional pressure regulators, the CS 215 thus contributing to the mobility of the entire Wenoll-System. This is also an important factor in mountain rescue and expedition medicinal uses. The weight of the complete CS 215 pressure regulator is just 670 g.

Periodic maintenance is at 5 year intervals, the date being documented on the test plaque.


The CS 215 pressure regulator is for administering medicinal oxygen in emergency medicine and for general O₂ therapy.

The multiple outlets allow two patients to be treated simultaneously.

Demand systems, breathing equipment or suction units can be attached to the medium pressure outlets.


■ Functional Components

- 1 High-pressure DIN 477 connector; 200 bar maximum
- 2 Cylinder pressure manometer
- 3 Medium pressure outlet, 5 bar, click connector
- 4 Constant flow outlet, 4-6 mm diameter, 0-20 l/min
- 5 Flush button with automatic reset (oxygen shower, 40 l/min)
- 6 Scale for oxygen flow, 0-20 l/min
- 7 Safety valve with protective cap
- 8 Test plaque with date of next safety check or maintenance
- 9 13 x 2.5 mm O-ring for high pressure connector

■ Flow Setting in Steps / Flush Function

Pressure regulator scale markings:

The flow setting is carried out in steps from 0-15 l/min. The fine settings up to 2.5 l/min are for applications in conjunction with a closed or recirculation system. The 10 and 15 litres/min. settings are used with an open system (O₂ mask with or without reservoir bellows). The oxygen flush button, with its automatic reset, is located at the centre of the adjustment knob. Operating this button applies a flushing flow of 40 litres/min. independent of the position of the knob.

■ Ordering Information

Pressure regulator CS 215, in steps 0-15 l/min, 200 bar, 2 medium pressure outlets (click connection) _____ Art.-Nr. 121 096

Medium pressure outlets also available with DIN (M12x1) and UNF (9/16"x18) - or combination.

The CS 215 pressure regulator is also available with British Standard winding (Art.-Nr. 121 097).

CO₂-Absorber with O₂ connecting tube

The CO₂-Absorber is the main item of the Wenoll re-circulation system. The application of the Wenoll re-circulation system provides multiple advantages:

- High effectiveness: 100% inspiratory O₂ concentration
- Ideal physiological conditions by breathing warm and humid oxygen
- Long therapy times (7.5 hours with a 2 litres O₂ cylinder, 200 bar)

Inspiration and expiration are automatically operated by valves, integrated in the absorber.

The connections for input/output as well as for the rebreathing bag are designed according to ISO standard. Furthermore they are colour marked. For the constantflow input the O₂ connecting tube is already assembled.

The absorber is filled with the absorbing medium Spherasorb. According to capacity of Spherasorb and the oxygen there were measured the overleaf application diagramms (depending on the used pressure regulator).


■ Functional components


- 1 Cartridge with sodalime Spherasorb 700g
- 2 Connector for tubing of breathing system (Yellow=expiration)
- 3 Connector for tubing of breathing system (green=inspiration)
- 4 Connector for rebreathing bag
- 5 O₂ connecting tube 1.6 m

■ O₂ dosage 1 patient with Wenoll re-circulation system and pressure regulator CS 215

■ O₂ dosage 1 patient with Wenoll re-circulation system and pressure regulator RS 220

■ Durability CO₂-Absorber

Due to hygiene factors the absorber has to be changed after operation (single use system). The storage period is 4 years, afterwards the absorber has to be replaced. There is a corresponding label with date of manufacture (LOT) and expiry date (red sandglass) on the absorber.

■ Dimensions and Weights CO₂-Absorber

Diameter: 145 mm
Weight, complete: 900 g

■ Ordering Information

CO₂-Absorber with O₂ connecting tube, 1.6 m _____ Art.-Nr. 121 050

New overflow valve as of 09/2022


The Wenoll-circulation system is overwhelmingly used for the spontaneous artificial respiration. In order to further facilitate the use of this system, the APL-valve 0-60 mbar used up to now was replaced by an overflow valve with a fixed setting at $P_{\text{max}} \sim 4$ mbar. This overflow valve is integrated into the expiration T-piece (yellow marking). The free ISO-fitting to this T-piece is protected by a white cover.

The APL-valve 0-60 mbar used up to now is hygienically packed and included in the standard scope of delivery of the Wenoll-cases WS 100/200/300.

For manual respiration, the white cover is removed and the APL-valve is fitted as usual.

Tubing system with overflow valve and cap for spontaneous respiration

Tubing system with overflow valve and APL valve for manual respiration

Wenoll Demand Module

The Wenoll demand module has been introduced as an option for the well-known Wenoll recirculation system.

The demand module consists of the demand valve with ISO connector and the oxygen connection tube. The demand valve itself was developed in cooperation with the Oceanic Company on the basis of the "SLIMLINE" octopus. The result was an oxygen demand valve matching medical pressure regulators while at the same time optimising patient breathing. The operating principle guarantees a concentration of 100% of the inhaled oxygen irrespective of how the patient is breathing, even if very shallowly.

A facemask or a mouthpiece may be used on the ISO 22/15 patient connector. Both of these components are contained as standard in the Wenoll emergency cases as delivered.

■ Usage

The Wenoll demand module is normally used for short-term O₂ therapy. This however assumes that the patient can breathe unassisted and is basically recommended for:

- Short therapy times
- Short transport times
- General circulatory problems
- Sufficiently large oxygen supply
- Two patients simultaneously

The patient connector is ISO 22/15, so that oxygen can be applied via mouthpiece, facemask or catheter. For reasons of hygiene and to improve breathing conditions, a HME* can be attached to the ISO connector of the O₂ demand valve. HMEs are used to physically warm and humidify the oxygen. The problems of hypothermia and dehydration are counteracted during normobar oxygen therapy; the patient breathes warm, humid oxygen.

* HME = Heat and Moisture Exchanger
(Art.-Nr. 121 052)

Important hint:

The Wenoll demand module is designed for the medical emergency at normal atmospheric pressure and may not be used under water.

■ Functional Components

- 1 O₂ demand valve
- 2 ISO 22/15 patient connector
- 3 O₂ flow, manual (oxygen shower)
- 4 O₂ medium pressure connector
- 5 O₂ connecting tube, 1.6 metres (click connector, DIN or UNF)


■ Breathing Action / Pressure Volume Diagram

The breathing action for inhalation and exhalation and the entire breathing action are shown in the accompanying diagram.

Breathing Action:

Out	0,38
In	0,31
Total	0,69
Joules/litre	

■ Retest

The Wenoll-demand module has to be retested every **5 years** at the manufacturer. There is a label with the date (month/year) of the next retest directly on the demand module.

■ Ordering information

- Wenoll demand module with O₂ connecting tube, 1.6 metres DIN _____ Art.-Nr. 121 070
- O₂ connecting tube, 1.6 metres with DIN thread _____ Art.-Nr. 121 071
- Wenoll demand module with O₂ connecting tube, 1.6 metres UNF _____ Art.-Nr. 121 072
- O₂ connecting tube, 1.6 metres with UNF thread _____ Art.-Nr. 121 073
- Wenoll demand module with O₂ connecting tube, 1.6 metres click connector _____ Art.-Nr. 121 074
- O₂ connecting tube, 1.6 metres with click connector thread _____ Art.-Nr. 121 075

Respiration Bag with Demand Module and O₂ Cylinder

The respiration bag MR-100 serves for the manual respiration of adults. The patient connection is executed in accordance with the ISO-standard 22/15. For the patients safety, an APL valve (60 mbar) is integrated. For an effective (100% O₂) respiration with an economic oxygen consumption, the respiration bag can be combined with the Wenoll-Demand module.

Respiration bag: Article no. 122 060

Bacterial filter / HME MEDISIZE HYGROVENT S

MEDISIZE HYGROVENT S is a breathing filter with the characteristics of an artificial nose which protects patients against bacterial and viral infections. MEDISIZE HYGROVENT S contains a hydrophobic filter which acts mechanically and electrostatically. Effectively warms and moistens the air to be inhaled, by means of the HME element and can in consultation with the doctor, replace the active humidifier. The connectors are according to ISO standards to assure the right assembly. CO₂-measurement via Luer Lock connector. MEDISIZE HYGROVENT S is equipped with a blue colored ring. In practice very easy to distinguish between the filter/HME and filter (blue=filter+HME). The ring also gives a better grip on the product during the assembly of the product into the system.

Warnings

- The HVS should be replaced immediately if contaminated by secretions or haemoptysis; this can result in significantly increased resistance.
- Do not use the HVS if the sterility indicator on the primary packaging is not yellow-brown in color.
- Do not use the HVS in combination with active humidifiers.
- Do not use the HVS in combination with explosive anaesthetic gases.
- Do not moisten the HVS prior to use - this must be done by the patient.
- Monitor the pressure drop across the HVS if it is being used in combination with medicine nebulisers.
- Maintenance of normal breathing pressure is important. This will prevent any condensation from being forced through the filter.
- Leakage around the endotracheal tube can cut down the effective humidification.

- The HVS should not be used if its packaging has been damaged.
- The HVS should only be used with patients who require humidification.

Utilization

- Check the sterility indicator on the primary packaging - it should be yellow-brown in color.
- The HVS should be fitted directly to the endotracheal tube or to a catheter mount and attached to the Y-connector of the breathing circuit.
- Before use check the breathing circuit for leakage.
- The HVS's maximum duration of use is 24 hours.

Technical Data

Bacterial retention	>99,99%
Virus retention	>99,99%
Mode of filtration	Elektrostatic / mechanical
Humidity of inspired air	33mg H ₂ O/l air at Vt 250ml (ISO 9360)
Pressure drop	<1,6 cm H ₂ O at 60 l/min
Weight	34g
Compressible volumen	55ml (ISO 9360)
Connectors	22M15F-22F/15M (ISO 5356-1)
Expiration in undamaged packaging	5 years if kept dry between 0°C and 40°C


Ordering Information

Bacterial filter / HME

Art.-Nr. 121 052

Adapter for O₂-cylinders DIN 477 / 200 bar

Type Nr. 1


suitable for following countries:

Egypt, United Kingdom,
Hongkong, India, Indonesia,
Irak, Brunei, Ireland, Jemen,
Jordania, Kuwait, Liberia, Libya,
Thailand, Malaysia, Malta,
Vietnam, Newsealand, Nigeria,
Pakistan, Saudi-Arabia,
Singapore, South Africa,
United Arab Emirates

Order information type 1:

Art.-Nr. 121 085

Type Nr. 2


suitable for following
countries:

Afghanistan, Argentina,
Brazil, Danmark, Finland,
Iran, Ethiopia, Iceland, Italia,
Chile, Ex-Jugoslawia,
Liechtenstein, Luxemburg,
Norway, Paraguay, Poland,
Albania, Sweden, Somalia,
Turkey, Czech Republic,
Hungary, Uruguay, Peru

Order information type 2:

Art.-Nr. 121 086

Type Nr. 3

suitable for following countries:

Algeria, Angola, Belgium,
Formosa, France, Greece,
Guinea, Libanon, Morocco,
Mosambique, Netherlands,
Portugal, Spain, Syra,
Tunesia

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- Service, Wartungen und Reparaturen in Zusammenarbeit mit EMS GmbH abzuwickeln.
EMS GmbH
Bevoegdheid tot nazicht en keuringen toegekend door EMS GmbH - Möhrendorf – Duitsland